ENGLISH WARNINGS

BEFORE USING THE OR CHARGING THE THERAGUN MINI, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL, ON THE CHARGER, AND ON THE Theragun mini. The Theragun mini is intended for commercial or home use. We strive to make the Theragun mini as safe as possible. This is an advanced mechanical device with electric components. If the Theragun mini and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not place your finger or any objects near the metal piston above the attachment while the Theragun mini is in use or immediately after use. The metal piston can become extremely hot. Do not use any attachments other than the ones provided by Therabody and only as instructed by Therabody. Do not use the Theragun mini on your head or near your genitals. Do not use the Theragun mini above the waist during pregnancy. Do not get hair caught in any moving part of the Theragun mini. Do not use the Theragun mini above your Adam's Apple or C4. Do not get the Theragun mini wet. Only clean the Theragun mini with a magic eraser sponge or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the Theragun mini. Never store the Theragun mini at temperatures above 40 C/104 F. Only charge the Theragun mini when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the Theragun mini or its charger to water. For long-term storage, the Theragun mini should be fully charged. This product is not intended to diagnose, treat or prevent any diseases. When using the Theragun mini, these basic precautions should always be followed including:

1. **USE ONLY AS INSTRUCTED.** Use the Theragun mini as described in this Theragun mini Instruction Manual. Use only Therabody recommended accessories and replacement parts. Do not carry out any maintenance other than what is shown in the Theragun mini Instruction Manual or advised by Therabody.
2. **NOT FOR CHILDREN.** The Theragun mini and charger are not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the Theragun mini to be used as a toy. Close attention is necessary when used by or near children. Children should be supervised to ensure that they do not play with the Theragun mini or charger.
3. **REDUCING RISK OF INJURY.** Charge the Theragun mini only with the Theragun mini charger. Other chargers may cause personal injury or damage. Do not connect the Theragun mini to a different power supply plug or car cigarette lighter. The internal battery will be permanently disabled or damaged.
4. **CHARGING LOCATIONS.** The Theragun mini should be charged indoors in a well ventilated, dry location. Do not charge the Theragun mini outdoors, in a bathroom or within 10 feet (3.1 metres) of a bathtub or pool.

Do not use the Theragun mini or charger on wet surfaces and do not expose the charger to moisture, rain or snow. Do not use the Theragun mini or charger in the presence of explosive atmospheres (gaseous fumes, dust or flammable materials). Sparks may be generated, possibly causing a fire.

5. **UNPLUGGING THE CHARGER.** Pull the plug, not the cable, to reduce the risk of damage to the electrical plug and cable. Never carry charger by its cable. Keep the cable from heat, oil and sharp edges. Do not use or charge with a damaged cable. Do not stretch the charger cable or place the cable under strain. Keep cable away from heated surfaces. Do not handle the charger, including the charger terminals or Theragun mini with wet hands.

6. **CHARGER CARE.** Unplug the charger when not in use. Unplug the charger from the socket when not in use for extended periods. Make sure the charger cable is located so that it will not be stepped on, tripped over, or otherwise subjected to damage or stress. Have a damaged charger replaced immediately.

7. **DO NOT OVERCHARGE.** Do not leave connect the Theragun mini to the charger for more than an hour after the Theragun mini has been fully charged. The Theragun mini includes a protection system to avoid the risk of overcharging, however, overcharging the Theragun mini may reduce its life over time. Only use with AC adapter: [HXD302-1501500CI/ZZ](#).

8. **DO NOT BLOCK MOTOR VENTS.** Do not put any obstructions into vents. Do not use if any vent is blocked. Keep free of dust, lint, hair and anything that may reduce airflow. Keep hair, loose clothing, and fingers from vents and moving parts.

9. **DO NOT BURN OR INCINERATE THE Theragun mini.** The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when battery is burned.

10. **DO NOT CRUSH, DROP OR DAMAGE THE Theragun mini CHARGER.** Do not use a charger that has received a sharp blow, been dropped, run over, or damaged in any way.

11. **BATTERY CHEMICALS CAUSE SERIOUS BURNS.** Never allow the internal battery to come into contact with the skin, eyes or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water and rinse with vinegar. If eyes are exposed to battery chemicals, immediately flush with water for 20 minutes and seek medical attention. Remove and dispose of contaminated clothing.

12. **DO NOT OPERATE UNDER BLANKET AND PILLOW.** Excessive heating can occur and cause fire, electric shock, or injury to persons.

13. **STORING THE Theragun mini AND CHARGER.** Store in a cool, dry place. Do not store the Theragun mini where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.

14. **DISPOSING OF BATTERIES.** The internal Theragun Lithium-ion batteries for the Theragun mini are more environmentally friendly than some other types of batteries (e.g. nickel-cadmium). Always dispose

of the Theragun mini according to government, regional and local regulations. Contact a recycling agency in your area for recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the Theragun mini's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.

15. **DO NOT DISASSEMBLE.** Disassembly or incorrect reassembly may result in the risk of electric shock, fire or exposure to battery chemicals. The warranty will be void if the Theragun mini or charger are disassembled or if any parts have been removed.

16. **SERVICE.** If the Theragun mini or charger is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use and contact Therabody by emailing info@therabody.com for repair or service. Do not attempt to repair or disassemble the Theragun mini which may result in an electric shock or fire.

17. **DO NOT USE WHILE BATHING OR IN THE SHOWER, BATHTUB OR SINK.** Do not place or store the Theragun mini where it can fall or be pulled into a bathtub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.

18. **THERMAL LIMITER.** The Theragun mini employs an automatic reset thermal limiter that shuts off the Theragun mini to prevent overheating and fire.

LIMITED WARRANTY

For full warranty information, please visit www.therabody.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address: Therabody - Warranty Attn: Customer Service 3003 W. Olympic Blvd., Suite 106-140 Los Angeles, CA 90006 Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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