CE FCC RoHS

(DE) (AT)

MASSAGEPISTOLE PRO 3

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7 (EN) INSTRUCTIONS

Congratulations on the purchase of this high-quality fitness device. This device has been examined and approved by recognized inspection and test facilities. These tests are carried out in the interests of the quality and safety of your training.

Please read this user manual carefully. The operating instructions are an integral part of the product. You should therefore keep these carefully, as well as the packaging, for later questions. If you pass the product on to third parties, please always include the operating instructions.

Please note that you have purchased a device for normal home use. This device is not intended for commercial use and is therefore not designed for loads such as those that occur in fitness studios.

8 SAFETY AND GENERAL INFORMATION

General Safety Information

- a) Please observe these operating instructions and keep them!
- b) Read and follow the instructions for use and all other information before using the device. Safe use of the device can only be achieved if it is used correctly. It is your responsibility to ensure that all users of the device are aware of the warnings and precautions.
- c) Be sure to consult your GP to identify and prevent any health risks before beginning any exercise program. Your doctor's advice is authoritative regarding your pulse rate, blood pressure or cholesterol level.
- d) Pay attention to your body signals. Incorrect or excessive training can endanger your health. If you experience pain, tightness in your chest, irregular heartbeat, extreme shortness of breath, dizziness or light-headedness at any time during your workout, please stop your workout immediately and consult your doctor.
- e) Keep children and animals away from the device. This device is intended for adult use only.
- f) Use the device only as specified. If you find a defective part during inspection, or if you hear unusual noises when using the device, stop using it until the problem is solved.
- g) The device is not suitable for therapeutic use.
- h) In the event of problems and/or a warranty claim, please send an email to SERVICE@LA-SPORTS.DE or use our contact form on the website <https://la-sports.de/kundenservice/> or contact us via WhatsApp on +49 2204 304616.
- i) Do not use the device continuously for more than one hour. Let the device rest for 30 minutes before using it again.
- j) Do not immerse the device in water. Keep it away from liquids or sources of heat. Keep the ventilation openings free from dust and dirt.
- k) Disconnect the charger from the battery of the device after charging or before use.

Specific Security Information

ATTENTION DANGER! To reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions:

- a) For use by adults only
- b) Use the product only on a dry clean surface of the body over clothing by gently pressing and moving across the body for approximately 60 seconds per region.
- c) Only use it on the soft tissues of the body without causing any pain or discomfort. Do not use on the head or any hard or bony area of the body.
- d) Only use the massage attachments that achieve the best possible effect.
- e) Bruising can occur regardless of control setting or pressure applied. Check the treated areas frequently and stop immediately at the first sign of pain or discomfort.
- f) Keep fingers, hair, or other body parts away from the barrel and back of the applicator as this may result in crushing.
- g) Do not place objects in the ventilation holes of the massage gun
- h) Do not immerse the massage gun in water or allow water to enter the massage gun ventilation holes
- i) Do not drop or abuse the massage gun.

- j) Only charge the device with the supplied charger.
- k) Check the device and battery carefully before each use.
- l) Do not tamper with or alter the massage gun in any way.
- m) Never leave the massage gun unattended
- n) Please do not use the massager without the express permission of a doctor if one or more of the following apply to you:
 - Pregnancy, diabetes with complications such as neuropathy or retinal damage, worn out pacemakers, recent surgery, epilepsy or herniated discs.
 - Spondylolisthesis, spondylolysis or spondylosis, recent joint replacements or IUDs, metal pins or plates, or concerns about your physical health.
 - Infirm persons should be accompanied by another adult when using a vibration device.
 - These contraindications do not mean that you cannot use vibration or exercise. However, we recommend that you consult a doctor first.

The effects of vibration exercise on certain medical disorders are currently being studied. This will likely shorten the list of contraindications shown above. Practical experience has shown that the integration of vibration exercises into a treatment plan is advisable in several cases. This must be done on the advice and under the supervision of a doctor, specialist, or physiotherapist.

Battery Safety Measures

- a) Do not short-circuit batteries.
- b) Do not heat batteries! Do not store your batteries near fire, stoves, or other heat sources. Do not place the batteries in direct sunlight. Do not use or store batteries in vehicles in hot weather.
- c) Do not weld or solder directly to batteries!
- d) Do not disassemble batteries!
- e) Do not deform batteries!
- f) Do not throw batteries into fire!
- g) Keep batteries out of the reach of children!

9 PACKAGING CONTENT

Before using your new device for the first time, please check that the individual parts are complete and intact. Please see figures A & B in the appendix. If, despite our quality controls, individual parts are missing, please contact our service center by email (SERVICE@LA-SPORTS.DE) using our contact form on the website: <https://la-sports.de/kundenservice/> or by WhatsApp at +49 2204 304616.

Rechargeable lithium-ion battery 13.5/26V
Mains voltage: 110-240V / 50/60Hz; 5A; 13.5/26VDC-1A

Extra light design for easy use (about 1.1 kg)

No.	Description
1	Powerful, high-torque brushless motor with Classic Quiet Smooth technology
2	Rechargeable lithium-ion battery (approx. 4 operating hours per charge)
3	Three (3) speed settings for optimal massage pressure
4	Vibration frequency: level 1: 30 Hz / level 2: 40 Hz / level 3: 50 Hz
5	On/Off switch (on the underside of the handle)
6	Four (4) interchangeable massage heads

10 USAGE INSTRUCTIONS

Charging

- (1) Before initial use, fully charge the battery for up to six (6) hours.
- (2) To charge, connect the appropriate plug of the supplied 24V adapter to the battery plug connector, and then plug the power adapter's power plug into the power outlet.
- (3) The LED strip on the battery lights will pulsate to indicate battery charge level and active charging.
- (4) Full battery charge is indicated when all LED indicators stop flashing and remain fully lit.
- (5) The battery can be charged at any time and at any battery level.
- (6) It is not recommended to fully discharge the battery down to the red LED.
- (7) Average run time is approximately two (2) hours depending on speed and pressure applied during use.

Usage/massage

- (1) Align and attach the battery to the base of the device. Make sure it snaps securely into place (if not pre-assembled).
- (2) With the device turned off, fully insert the desired massage attachment. Pull it out (also with the device turned off) by pressing remove while pulling it out.
- (3) Turn on the device by turning the power switch on the bottom of the battery to ON. Turn off the device by turning the power switch to OFF.
- (4) With the power on, press the button once to turn on the device at level one, twice to turn on the device at level two, and three times to turn on the device at level three. Pressing the power button again will turn off the device. The blue LED lights indicate the corresponding power level.
- (5) Massage the desired body part at the desired speed, applying more or less pressure as needed.
- (6) When changing the massage head, you may apply some vegetable oil to the inner wall of the massage head or the outer diameter of the massage head rubber ring to increase the degree of lubrication of the massage head and the mounting hole and make it easier to take out. It is recommended to store the massage head separately to prevent sticking to the massage gun mounting hole.

11 MAINTENANCE AND CARE

General information

Using your device will naturally cause wear and tear. This requires regular care and maintenance. Under normal conditions and use, your device is maintenance-free. However, it is the user's responsibility to periodically inspect these parts for wear and tear that will affect performance.

Note that the performance of your device can be significantly reduced by the influence of environmental conditions. Do not expose your device to excessive heat or extreme cold for long periods of time as performance may be affected.

Any modification to the original parts, materials, or construction of the device can significantly reduce performance. The same applies in the event of heavy soiling or misuse.

Do not use aggressive cleaning agents for cleaning, but remove dirt and dust with a soft, slightly damp cloth.

Disposal Instructions

- (1) This product is shipped in packaging to protect against shipping damage. The packaging materials are made from raw materials that can be recycled. Please dispose of these materials separately.
- (2) If you wish to part with this product, dispose of it according to current regulations. Information can be obtained from the municipal office.
- (3) Batteries do not belong in household waste! Dispose of batteries in an environmentally friendly manner and dispose of them properly in the interest of the environment. Your local administration will inform you about collection points and opening times.

12 WARRANTY CONDITIONS

L.A. Sports GmbH & Co. KG grants you a two (2) year guarantee for this product from the date of purchase in accordance with our warranty provisions and the applicable legal regulations. This warranty applies to first-time retail purchasers only and is not transferable.

Despite the extensive quality controls and the use of high-quality materials, defects in this product can occur in exceptional cases. Within the guarantee period, we will rectify any defects that can be proven to be due to material or manufacturing defects free of charge. Further or other claims, in particular for compensation for damage not caused to the product itself, are excluded.

The warranty does not cover damage caused by natural wear and tear, commercial use, improper use, incorrect operation, lack of care, incorrect installation and storage, non-compliance with the operating instructions, as well as force majeure or other external influences. The guarantee is void if repairs or interventions are carried out by persons who are not authorized to do so.

Warranty claims must be asserted within the warranty period. The fully completed warranty certificate and proof of purchase must be enclosed with the device. The device must always be sent to the service center in the original packaging.

To order spare parts, please contact our service center using the following contact options, thank you:

- by email to SERVICE@LA-SPORTS.DE
- via our contact form at <https://LA-SPORTS.DE/KUNDENSERVICE>
- via WhatsApp message to 02204 - 304616 (scanning the adjacent QR code with your smartphone will lead you directly to a WhatsApp message to our service center)

For other warranty claims or faults, please contact your sales outlet.

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Név: L.A. Sports GmbH & Co. KG